

SONIMAGE P3

Personal, Portable, Point-of-Care Ultrasound

Konica Minolta Medical Imaging introduces the new SONIMAGE P3, a personal, portable, point-of-care ultrasound system that revolutionizes how you see your patients.



No Need to Compromise. Ultras

PERSONAL

SONIMAGE P3 strengthens clinical confidence, aids in quick diagnosis, and enhances the patient and healthcare professional experience. The ultrasound technology is easy to use in various clinical settings.

PORTABLE


SONIMAGE P3 is small and lightweight, enabling imaging at the point-of-care. It powers up instantly and is ready to use when and wherever you need it. Use as a stand-alone system or simply plug the handheld display unit or ultrasound probe directly into your Windows based PC, laptop or tablet running SONIMAGE VIEW software – standard with the SONIMAGE system. Ideal for occasions when you need a larger display.

POINT-OF-CARE

SONIMAGE P3 is designed as a handheld system for emerging point-of-care markets, with a broad range of applications for both human and veterinarian markets. Ideal for any point-of-care physician or healthcare professional.

BROAD RANGE OF APPLICATIONS

- FAST and eFAST
- Lung imaging pleural fluid, thoracic and pleural motion
- Post-void residual bladder volume (automated measurements)
- Basic pregnancy assessments confirmation scans, FHR, fetal position
- Basic internal organ imaging gross gall bladder pathology
- Abdominal aortic aneurysm (AAA) Screening
- Localization of free fluid in the abdomen and lungs
- Core muscle assessment Transverse Abdominus, pelvic floor



Abdomen

"In my experience, B line, C line and pleural fluid resolution is superb."

Dr. Pat Allan

ound Designed for You.

INCREASE EFFICIENCY. IMPROVE WORKFLOW.

SONIMAGE P3's quick diagnosis during the examination helps reduce unnecessary follow-ups and can improve referral accuracy. Immediate ultrasound helps speed diagnosis, reduce patient wait-times, and improve work flow by providing additional information __ for validation and identification of patients who require further examination.



Innovative design makes this real-time system an ideal companion for any point-of-care physician or healthcare professional. The convenience of B-mode, M-mode, Doppler and simple to access presets, makes SONIMAGE P3 simply indispensable.





Personal, Portable, Point-of-Care Ultrasound

A revolution improving the Quality-of-Care at the Point-of-Care

"Dozens of uncomplicated thoracentesis have been done. Unlike its competitors, it is very easy to transfer images from the SONIMAGE P3 using either direct USB or "jump" drive storage. These images were often accompanied by digital and audio annotations that were well-suited to medical education purposes."

Dr. Pat Allan, Pulmonologist Dayton Respiratory Center, Dayton, OH

FEATURES	
Weight	• 392 grams (13.8 oz)
Dimensions	 Display unit: 11.5cm H x 6.8cm W x 2.0cm D Display Screen: 6.85cm Hi-res colour LCD Probe: 14cm L x ~6cm W x ~3.3cm D Connecting cable: 85.8cm L
Transducer	• S3-5 Sector Annular (3-5 MHz) – B/M Mode + PW Doppler
User interface	 Touch screen and scroll wheel Onscreen alphanumeric keyboard Text and Voice Recording Optional tablet display PictBridge compatible printing
Measurement tools	 Auto Volume Measurements Linear distance calipers Ellipse & polygon area/circumference Beats-per-minute (BPM)
Image storage	4 GB microSD data card (included)Urgent Mode internal memoryUSB connection for downloads
Software	 SONIMAGE VIEW PC/tablet software Backup images/data Search and manage database Automatic software updates
Battery	 75% charge in 2 hours 100% charge in 3 hours Lasts up to 30 hours (based on four 1-minute scans/hour) Over 120 minutes of continuous scanning



KONICA MINOLTA MEDICAL IMAGING 411 NEWARK POMPTON TURNPIKE WAYNE, NJ 07470

TEL: (973) 633-1500 FAX: (973) 523-7408 WEBSITE: medical.konicaminolta.us